

426

Advanced sample preparation

Mr Pawel Nowakowski

EA Fischione Instruments Inc, Pittsburgh, USA

Abstract

Fischione Instruments' Model 1062 TrionMill is an excellent tool for creating the sample surface characteristics needed for SEM imaging and analyses.

The TrionMill's three ion sources provide large-scale milling of planar and cross-section samples. High-energy operation allows rapid milling of large areas (up to 50 mm); low-energy operation allows gentle sample polishing.



FISCHLINE
LABORATORY

Mo